

Features

- Digital timer adjustable upto 99 min.
- Replaceable bulbs
- Stainless steel chamber
- Compact Size

Technical Details


- Dimension in mm: 250 X 150 X 200
- Weight: 1.125kg
- Mains Supply: 230 Volts, 50/60Hz
- Power Consumption: 22 Watts
- Wavelength: 450nm



Available at:

Delta

Blu-lux



Delta

Manufactured by:

Delta Labs

323 SIDCO Women Industrial Estate, Vellanur P.O


Thirumullaivoyal, Chennai - 600 062. TN. INDIA

T: +91 98849 22713

E: dentalabsindia@gmail.com

W: deltadentallabs.com

Blu-lux is a light curing unit used for curing light curing trays for the fabrication of custom impression trays. It is also used for post curing of 3D printed resins and light curing pattern resins.



Special tray for complete dentures with Intra Oral tracer attached



Custom tray for implant impressions



Poly Tray

Poly Tray is light curing tray material for the fabrication of individual trays or special trays for crown and bridge, implants and complete dentures.

Features

- Easy handling and adaptation
- Precise fitting
- Economical

Article No. 10074 Qty : 50 Pcs

